

Process Set-Up

1. Turn **Thermolyne** digital hotplate **ON** (Rm 201 K) and set temperature to 110°C. Once set point is reached, verify temperature measures 90°C +/- 2C° using **Omega HH508** digital thermometer/probe.
2. Turn on **DI Water** in service bay (Rm 201 D).
3. Turn **Quintel Q-4000** mask aligner vacuum **ON**. The vacuum switch is located on the back wall to the left of mask aligner. If the **HTG** mask aligner in room 201 H is used, verify vacuum is at least 20 in. Hg as indicated on the vacuum gauge located behind the **HTG** mask aligner.
4. Turn mask aligner power **ON** and verify UV lamp is operational. Measure the UV intensity using the **Quintel Vari-Wave II** meter prior to exposure. Calculate exposure time using the following:
$$\text{Exposure Time (seconds)} = \text{Exposure (mJ/cm}^2\text{)} / \text{Intensity (mW/cm}^2\text{)}$$
5. Turn exhaust to spinner bowl (room 201 K) **ON** by connecting vacuum plug to power outlet on front left bottom of exhaust hood.
6. Locate **JSR NFR-014R** photoresist, **PD 523AD** developer solution, and **Fisher** disposable polyethylene pipettes.

Process

Spin Coat (Headway Spinner - Rm 201 K) Test spin for speed before applying photoresist and adjust if necessary. Apply **JSR NFR-014R** resist using a disposable pipette onto wafer. Spin @ 3K RPM for 60 seconds.

Soft Bake on hotplate @ 90°C (set point = 110°C) for 90 seconds.

Expose (HTG or Quintel) Insert proper mask into exposure system, chrome side down. Expose at 80 mJ/cm². Set exposure time based on intensity calculation (See setup procedure for calculation). Align to previous pattern if necessary.

Post Exposure Bake on hotplate @ 90°C (set point = 110°C) for 90 seconds

Develop Immerse wafer(s) into developer solution **PD 523AD**, full strength for 60 seconds with constant agitation (one swish per second). Dispose of developer solution once develop is completed. Rinse in **DI H₂O** for 5 minutes and transfer to cascade rinse tank in room 201 J. Rinse wafers until resistivity is >15 MW. Transfer wafers into **STI Semitool** spin dryer and press **START**. **Note:** Use 1000 mL **PD 523AD** developer for 1-3 wafers, 2000 mL **PD 523AD** developer for 4-6 wafers, etc.

Inspect pattern. Proceed to next step if OK. Strip photoresist and return to spin coat if unacceptable.

Hard Bake Bake in convection oven for 20 minutes @ 90°C. **Note:** this is optional if film is thin (1-2Å).

Etch using the appropriate etch for the material to be etched.

Inspect for completeness of etch. Etch for additional time if necessary.

Strip Photoresist using either standard single use **Piranha** solution or **CARC Resist Remover** solution.