

**Typical Procedure**

- Preheat hot plates to 60°C and 100°C.  
**Note:** See DataPlate® Hotplate manual.
- Clean sample in acetone, methanol, isopropanol solution at 60°C for 5 min.
- During this time,
  - Prepare two small baths (~25 ml each) of 1:4 [400K developer:DI H<sub>2</sub>O].
  - Start mask aligner. See Karl Suss Mask Aligner manual.
  - Start spincoater. See SCS SpinCoater manual.
  - Program appropriate spin parameters as recommended at right.
- Remove solvent from hotplate; change hotplate temperature to 150°C.
- Dry sample with compressed N<sub>2</sub>; clean in diluted developer by dipping for 40 sec in first bath, 20 sec in second bath.  
**Note:** Second bath will be cleaner and may be used later. Dispose of first bath. Prepare small bath (~25ml) for later use, and set both aside.
- Rinse sample in acetone, then methanol, then isopropanol.
- Dry with compressed N<sub>2</sub>.
- Place sample on spincoater chuck. (A test spin is recommended.)
- Cover sample completely with LOR (10B or 30B).
- Press Start button and allow spin recipe to complete.  
**Note:** If resist appears on backside of sample, drop 2-3 drops of EBR remover on clean room towel and carefully glide backside of sample across wetted area. Take care not to expose topside to remover.
- Bake on 150°C hotplate for 5 min; change hotplate temperature to 110°C.
- Cool sample on cold surface for 10 sec.
- Spin-test sample on spincoater.
- Cover sample completely with Clariant AZ 1518 photoresist .
- Press Start button, and allow entire spin recipe to complete.  
**Note:** If resist appears on backside, clean with EBR remover
- Soft bake on 100°C hotplate for 45 sec.
- Cool sample on cold surface for 10 sec.
- Select chrome or FeO mask and rinse pattern-side of mask with acetone.
- Dry mask thoroughly with compressed N<sub>2</sub>.
- Develop sample by dipping in 1:4 developer for 40 sec, then 20 sec in a second bath.
- Rinse with DI water and dry with compressed N<sub>2</sub>.
- Inspect photoresist pattern quality using a microscope.
- If pattern quality is good, hard bake on 110°C hot plate for 5 minutes.  
**Note:** Further develop in 10 sec intervals if needed. Otherwise, strip PR and LOR thoroughly with acetone and begin again.
- Remove sample from hotplate and cool on cold surface before storing in carrier.
- Turn off hotplates, aligner, spincoater and microscope.

**Clariant AZ 1518**

RPM1 = 2000  
RAMP1 = 2  
TIME1 = 1  
RPM2 = 3000  
RAMP2 = 1  
TIME2 = 1  
RPM3 = 4000  
RAMP3 = 1  
TIME3 = 30

**LOR-10B**

RPM1 = 2000  
RAMP1 = 2  
TIME1 = 1  
RPM2 = 3000  
RAMP2 = 1  
TIME2 = 1  
RPM3 = 4000  
RAMP3 = 1  
TIME3 = 30

**LOR-30B**

RPM1 = 500  
RAMP1 = 2  
TIME1 = 1  
RPM2 = 900  
RAMP2 = 5  
TIME2 = 1  
RPM3 = 1000 – 1500  
RAMP3 = 1  
TIME3 = 30 – 60

**Clean-Up**

- While photoresist is hard-baking, remove glass mask from mask chuck and clean patterned side with acetone to remove any residual PR. Dry with compressed N<sub>2</sub> gas, and return to protective case.
- Clean up wet bench by disposing of solvents and developer.  
**Note:** Never leave open solvent waste bottle in sink. This is a violation of chemical safety procedure.
- Rinse glassware and used containers with DI H<sub>2</sub>O, return upside down on drying rack.
- Dry bench top.